



NGS TP Index Kit for Illumina (Index Primers Set II)

转座酶法二代测序多样本引物试剂盒 II (Illumina)

目录号：CW2959S (384 rxns)

保存条件：-20℃保存, 干冰运输。

产品内容

Component	CW2959S
	384 rxns
Index N502-N522 Primers for Illumina	16×24 μl
Index N701-N729 Primers for Illumina	24×16 μl

产品简介

本试剂盒是转座酶法二代测序快速DNA建库试剂盒专用的配套试剂盒，专为Illumina平台建库设计，其中包含16种N5端引物和24种N7端引物，可用于制备384种不同的双端Index文库。试剂盒中提供的所有试剂都经过严格的质量控制和功能验证，最大程度上保证了文库构建的稳定性和重复性。制备的文库可用于HiSeq X-10/4000/2500/2000和MiSeq等Illumina平台测序。

自备仪器、试剂和耗材

1. 磁力架：建议使用DynaMag™-2 (Cat.No. 12321D)。
2. DNA纯化回收试剂盒：建议使用康为磁珠法DNA纯化回收试剂盒 (Cat.No. CW2508)。
3. DNA建库试剂盒：建议使用康为世纪转座酶法二代测序快速DNA建库试剂盒 (Cat.No. CW2845/CW2846/CW2847)。
4. 无水乙醇。
5. 反应管：建议使用低吸附的PCR管与1.5 ml离心管；
枪头：建议使用高质量过滤枪头防止试剂盒、文库样本污染。

实验前准备及重要注意事项

开盖前请短暂离心，使液体收集到管底，以免不同引物间交叉污染。

操作步骤

康为世纪转座酶法二代测序多样本引物试剂盒使用方法请按照康为世纪转座酶法二代测序快速DNA建库试剂盒（Cat .No.CW2845/CW2846/CW2847）protocol进行。

Index N502-N522 Primers for Illumina

	Index Primers for Illumina	Index
N502	5'-AATGATACGGCGACCACCGAGATCTA CAC CTCTCTAT TCGTCCGGCAGCGTC-3'	CTCTCTAT
N503	5'-AATGATACGGCGACCACCGAGATCTA CA TATCCTCT TCGTCCGGCAGCGTC-3'	TATCCTCT
N505	5'-AATGATACGGCGACCACCGAGATCTA CAC GTAAGGAG TCGTCCGGCAGCGTC-3'	GTAAGGAG
N506	5'-AATGATACGGCGACCACCGAGATCTA CAC ACTGCATA TCGTCCGGCAGCGTC-3'	ACTGCATA
N507	5'-AATGATACGGCGACCACCGAGATCTA CAC AAGGAGTA TCGTCCGGCAGCGTC-3'	AAGGAGTA
N508	5'-AATGATACGGCGACCACCGAGATCTA CAC CTAAGCCT TCGTCCGGCAGCGTC-3'	CTAAGCCT
N510	5'-AATGATACGGCGACCACCGAGATCTA CAC CGTCTAAT TCGTCCGGCAGCGTC-4'	CGTCTAAT
N511	5'-AATGATACGGCGACCACCGAGATCTA CAC TCTCTCCG TCGTCCGGCAGCGTC-5'	TCTCTCCG

	Index Primers for Illumina	Index
N513	5'-AATGATACGGCGACCACCGAGATCTA CACT CGACTAG TCGTCGGCAGCGTC-6'	TCGACTAG
N515	5'-AATGATACGGCGACCACCGAGATCTA CACT TTCTAGCT TCGTCGGCAGCGTC-7'	TTCTAGCT
N516	5'-AATGATACGGCGACCACCGAGATCTA CAC CCTAGAGT TCGTCGGCAGCGTC-8'	CCTAGAGT
N517	5'-AATGATACGGCGACCACCGAGATCTA CAC GCGTAAGA TCGTCGGCAGCGTC-3'	GCGTAAGA
N518	5'-AATGATACGGCGACCACCGAGATCTA CAC CTATTAAG TCGTCGGCAGCGTC-4'	CTATTAAG
N520	5'-AATGATACGGCGACCACCGAGATCTA CAC AAGGCTAT TCGTCGGCAGCGTC-5'	AAGGCTAT
N521	5'-AATGATACGGCGACCACCGAGATCTA CAC GAGCCTTA TCGTCGGCAGCGTC-6'	GAGCCTTA
N522	5'-AATGATACGGCGACCACCGAGATCTA CACT TTATGCGA TCGTCGGCAGCGTC-7'	TTATGCGA

Index N701-N729 Primers for Illumina

	Index Primers for Illumina	Index
N701	5'-CAAGCAGAAGACGGCATAACGAGAT TCGCCTTA GTCTCGTGGGCTCGG-3'	TAAGGCGA
N702	5'-CAAGCAGAAGACGGCATAACGAGAT CTAGTACG GTCTCGTGGGCTCGG-3'	CGTACTAG
N703	5'-CAAGCAGAAGACGGCATAACGAGAT TTCTGCCT GTCTCGTGGGCTCGG-3'	AGGCAGAA
N704	5'-CAAGCAGAAGACGGCATAACGAGAT GCTCAGGA GTCTCGTGGGCTCGG-3'	TCCTGAGC
N705	5'-CAAGCAGAAGACGGCATAACGAGAT AGGAGTCC GTCTCGTGGGCTCGG-3'	GGACTCCT
N706	5'-CAAGCAGAAGACGGCATAACGAGAT CATGCCTA GTCTCGTGGGCTCGG-3'	TAGGCATG
N707	5'-CAAGCAGAAGACGGCATAACGAGAT GTAGAGAG GTCTCGTGGGCTCGG-3'	CTCTCTAC
N710	5'-CAAGCAGAAGACGGCATAACGAGAT CAGCCTCG GTCTCGTGGGCTCGG-3'	CGAGGCTG

	Index Primers for Illumina	Index
N711	5'-CAAGCAGAAGACGGCATAACGAGAT TGCCTCTT GTCTCGTGGGCTCGG-3'	AAGAGGCA
N712	5'-CAAGCAGAAGACGGCATAACGAGAT TCCTCTAC GTCTCGTGGGCTCGG-3'	GTAGAGGA
N714	5'-CAAGCAGAAGACGGCATAACGAGAT TCATGAGC GTCTCGTGGGCTCGG-3'	GCTCATGA
N715	5'-CAAGCAGAAGACGGCATAACGAGAT CCTGAGAT GTCTCGTGGGCTCGG-3'	ATCTCAGG
N716	5'-CAAGCAGAAGACGGCATAACGAGAT TAGCGAGT GTCTCGTGGGCTCGG-3'	ACTCGCTA
N718	5'-CAAGCAGAAGACGGCATAACGAGAT GTAGCTCC GTCTCGTGGGCTCGG-3'	GGAGCTAC
N719	5'-CAAGCAGAAGACGGCATAACGAGAT TACTACGC GTCTCGTGGGCTCGG-3'	GCGTAGTA
N720	5'-CAAGCAGAAGACGGCATAACGAGAT AGGCTCCG GTCTCGTGGGCTCGG-3'	CGGAGCCT
N721	5'-CAAGCAGAAGACGGCATAACGAGAT GCAGCGTA GTCTCGTGGGCTCGG-3'	TACGCTGC
N722	5'-CAAGCAGAAGACGGCATAACGAGAT CTGCGCAT GTCTCGTGGGCTCGG-3'	ATGCGCAG
N723	5'-CAAGCAGAAGACGGCATAACGAGAT GAGCGCTA GTCTCGTGGGCTCGG-3'	TAGCGCTC
N724	5'-CAAGCAGAAGACGGCATAACGAGAT CGCTCAGT GTCTCGTGGGCTCGG-3'	ACTGAGCG
N726	5'-CAAGCAGAAGACGGCATAACGAGAT GTCTTAGG GTCTCGTGGGCTCGG-3'	CCTAAGAC
N727	5'-CAAGCAGAAGACGGCATAACGAGAT ACTGATCG GTCTCGTGGGCTCGG-3'	CGATCAGT
N728	5'-CAAGCAGAAGACGGCATAACGAGAT TAGCTGCA GTCTCGTGGGCTCGG-3'	TGCAGCTA
N729	5'-CAAGCAGAAGACGGCATAACGAGAT GAGCTCGA GTCTCGTGGGCTCGG-3'	TCGACGTC